


## Dräger Breathing Bags Consumables and Accessories

Dräger Breathing Bags – It is your choice! More comfort and safety in anaesthesia: We offer breathing bags as a separate bag, as VentStar® Bag Set and as Infinity® ID BagSystem (with RFID-technology). The latex-free breathing bags are available in four sizes to fit any patient.

- 4 sizes to fit any patient
- Single patient use
- Handling by customer
- Bio-compatibility testing



## Benefits

---

### Four sizes to fit any patient

The latex-free separate breathing bags are available in four sizes to fit any patient. We offer breathing bags for:

- adult ( 2 and 3 L)
- pediatric ( 1 L)
- neonatal (0.5 L)

---

### Hygienic and practical: breathing bag for single use

Dräger breathing bags meet the hygiene requirements and help prevent cross-contamination. The single use not only ensures safety, but also contributes to process efficiency by allowing caregivers to allocate more time to other aspects of care.

---

### Safe handling by Customer

The breathing bags are made of a durable, compliant membrane that withstands demanding workflow. Dual-pleated sides and surface texture provide maximum grip.

---

### Bio-compatibility testing

The breathing bags are not made with natural rubber latex. The bags had been tested for their bio-compatibility. All Components are free of PVC and do not contain phthalate (DEHP) still Bisphenol A (BPA).

## Technical Data

### Sizes

Neo	Pediatric	Adult
0.5 L	1.0 L	2.0 L, 3.0 L
Material		Breathing bag hose EVA
		Connectors PP
		Hose holder TPE
		Breathing bag chloroprene, PP

## Ordering Information

### Breathing Bag Portfolio – Adult

VentStar® Breathing Bag 3 L	Packing Unit: 25	Bag only	MP02713
VentStar® Breathing Bag 2 L	Packing Unit: 25	Bag only	MP02712
VentStar® Breathing Bag Set 2 L	Packing Unit: 25	Bag 2 Liter; Hose 0.8 meter	MP00222
VentStar® Breathing Bag Set 2 L	Packing Unit: 25	Bag 2 Liter; Hose 1.2 meter	MP00380
VentStar® Breathing Bag Set 2 L	Packing Unit: 25	Bag 2 Liter; Hose 1.5 meter	MP00381
Infinity® ID Breathing Bag Set 2 L	Packing Unit: 25	Bag 2 Liter; Hose 1.2 meter, RFID Tack	MP00341
Infinity® ID Breathing Bag Set 2 L	Packing Unit: 25	Bag 2 Liter; Hose 1.5 meter, RFID Tack	MP01330

### Breathing Bag Portfolio – Paediatric

VentStar® Breathing Bag 1 L	Packing Unit: 25	Bag only	MP02711
VentStar® Breathing Bag Set (P)	Packing Unit: 25	Bag 1 Liter; Hose 1.1 meter	MP00384
VentStar® ID Breathing Bag Set (P)	Packing Unit: 25	Bag 1 Liter; Hose 1.1 meter, RFID Tack	MP01343

### Breathing Bag Portfolio – Neonatal

VentStar® Breathing Bag 0,5 L	Packing Unit: 25	Bag only	MP02710
VentStar® Breathing Bag Set (N)	Packing Unit: 25	Bag 0.5 Liter; Hose 1.1 meter	MP00383
VentStar® ID Breathing Bag Set (N)	Packing Unit: 25	Bag 0.5 Liter; Hose 1.1 meter, RFID Tack	MP01329

Not all products are worldwide available.

## Notes

### CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
[www.draeger.com](http://www.draeger.com)

### Manufacturer:

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)



### REGION DACH

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION EUROPE

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.  
25 International Business Park  
#04-20/21 German Centre  
Singapore 609916  
Tel +65 6308 9400  
Fax +65 6308 9401  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### REGION NORTH AMERICA

Draeger, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 800 4DRAGER  
(+1 800 437 2437)  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

### REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial  
S. de R.L.  
Complejo Business Park,  
V tower, 10th floor  
Panama City  
Tel +507 377 9100  
Fax +507 377 9130  
[contactcsa@draeger.com](mailto:contactcsa@draeger.com)