

Carrying System

This ergonomically designed Carrying System provides an integrated solution for carrying and transporting an Oxylog® 3000, Oxylog® 2000 plus or an Oxylog® 3000 plus and an oxygen cylinder with a minimum of effort. The frame is small and light with rounded corners to prevent any discomfort in carrying.



MT-1827-2007

When a patient is being transported, the complete Carrying System can simply be hung onto a variety of beds and stretchers or onto standard rails.

The bag can be used for storing all kinds of accessories: the patient ventilation tube can be stored in the front bag, even when already connected to the Oxylog 3000, Oxylog 2000 plus or to the Oxylog 3000 plus.

To mount the Carrying System in the All-Round Wall Holder, all you have to do is push one handle. The Wall Holder also

has an integrated automatic power connection so you no longer need to manually connect a DC/DC converter or AC/DC power supply.

With the Automatic Gas Source Switch, automatic switch-over from central gas supply to cylinder operation and vice versa becomes possible, eliminating the need to reconnect gas hoses. During switch-over ventilation continues as normal, thus offering smooth therapy.



MT-4253-2007

Complete ventilation solution in just one hand

SPECIFICATIONS

Size (width × height × depth)	515 × 400 × 260 mm / 20.2 × 15.7 × 10.2 inch (standard cylinder bracket)
Weight including Oxylog 3000, Oxylog 2000 plus or Oxylog 3000 plus 2L GFC cylinder, pressure reducer, oxygen connection tubes.	approx. 14,5 kg / 32 lbs
Cylinder sizes	
Diameter	90 mm to 120 mm / 3.5 to 4.7 inch
Length (cylinder bracket: standard)	300 mm to 350 mm / 11.8 to 13.8 inch (excl. valve)
Length (cylinder bracket: long)	300 mm to 500 mm / 11.8 to 19.7 inch (excl. valve)

ORDER LIST

Carrying System	2M 86 975
For carrying, transporting and operating Oxylog 3000, Oxylog 2000 plus or an Oxylog 3000 plus when no external power supply or central gas supply is available. Including carrying frame, accessories bag, dapter for quick power connection and cylinder bracket available in two lengths.	
To choose	
Cylinder bracket standard (1 – 2.5 L) For cylinders up to 120 mm / 4.7 inch diameter and 350 mm /13.8 inch length.	57 04 449
Cylinder bracket long (1 – 3 L) For cylinders up to 120 mm / 4.7 inch diameter and 500 mm / 19.7 inch length.	57 04 589
Oxygen supply	
Gas Supply System	57 045 00
To connect the Oxylog with an oxygen cylinder and/or a central gas supply. The central gas supply is automatically selected by the optional Automatic Gas Source Switch when the oxygen cylinder and central gas supply are simultaneously connected. Including oxygen connection tubes and optional pressure reducer.	
Optional accessories	
All-Round Wall Holder	57 04 216
Quick Power Connector	57 04 217

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany

As of August 2015:

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA.

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com