



Reliability goes XL

DRÄGER
FABIUS® PLUS XL

“Performance you can **count on**”



RELIABILITY IS KEY

There are many advanced anesthesia devices on the market today. And even the most advanced technology has to measure up on a day-to-day basis. The reliability of medical technology is a key factor for successful therapy as well as for economic success in an increasingly dynamic and competitive health care market.

SIMPLE IS BETTER

For safety and efficiency reasons medical technology should be simple to operate. Often, your personnel has little time to become familiar with all the functions and capabilities of a medical device so that they can react quickly and correctly in emergency situations. Caregivers need information that is well organized, highly visible and easy to read.

FLEXIBILITY IS A MUST

Hospitals should also be able to count on built-in expandability and competent, dedicated manufacturer support which begins from the first day and continues until the device is taken out of service. This helps caregivers to provide the best possible treatment over time.



THE FABIUS PLUS XL

Designed to provide quality therapy and sophisticated monitoring capabilities in a compact yet versatile package, the Fabius plus XL incorporates proven German engineering. With a range of expansion options to choose from, it lets you enjoy the advantages of Dräger quality today while maintaining flexibility for the future.

POWERFUL YET AFFORDABLE

In addition to the high-quality ventilation capabilities made possible by its highly accurate, electronically driven piston ventilator, the Fabius plus XL includes a number of features which help you maximize your therapeutic potential and improve patient safety.

A CLEAR VIEW

The Fabius plus XL organizes and displays information clearly and logically on a 10.4 inch color display, providing an excellent overview from a wide range of viewing angles.

FLEXIBLE SCALABILITY

The Dräger Fabius plus XL lets you flexibly expand its capabilities by mounting additional monitors, gas cylinders and accessories such as infusion pumps.

IMPROVED ERGONOMICS

Carefully positioned mounting rails let you mount additional equipment within easy reach. The optional higher trolley lets you work comfortably in a standing position. An enhanced cable management combined with integrated power sockets mean less cable clutter.

THE SPACE YOU NEED

Three large drawers let you store consumables where they are needed – close by. The top storage drawer can be converted to install your choice of Scio gas measurement module.

A CLEAR VIEW

- 10.4" color display
- Precise, mechanical gas flow tubes
- Visible ventilator movement
- Front visible pressure gauges for mounted gas cylinders

OPTIONS FOR FLEXIBLE SCALABILITY

- Scio or Vamos gas measurement integration
- Various mounting options of patient monitors
- Infusion pump mounts
- Integrated, medical grade power sockets

IMPROVED ERGONOMICS

- Enhanced cable and hose management
- Leak tight, compact breathing system with optional heater
- Vapor mount for 2 vapors
- Optional higher trolley
- Available as trolley, ceiling or wall-mounted version

THE SPACE YOU NEED

- Large drawers for consumables
- Improved storage space on top plate
- Optional basket for even more storage space



D-47654-2012

Ready for tomorrow



VENTILATE WITH ENHANCED CONFIDENCE

The Fabius plus XL is equipped with our powerful and precise E-vent piston ventilator, which provides high tidal volume accuracy and requires no drive gas for operation. This important advantage lets you continue to use it temporarily should the central gas supply fail. It comes complete with advanced features such as fresh gas decoupling, which maintains constant tidal volumes regardless of fresh gas flow changes. Delivering ICU-like ventilation performance and supporting all major ventilation modes, the Fabius plus XL is ready for a broad range of ventilation applications and acuity levels.

ECONOMICAL OPERATION

With its optional breathing circuit heater, the Fabius plus XL is an excellent choice for economical low-flow ventilation. This ventilation strategy reduces the consumption of costly anesthetic agents. The highly reliable Fabius plus XL helps maximize device availability, letting you treat more patients at less cost.

EASE OF USE

Thanks to its simple operating concept and intuitive menu architecture, the Fabius plus XL is easy to use and understand. Training efforts can be kept to a minimum, which can help reduce user errors.

Build your own Fabius plus XL



WORKSTATION WITH INFINITY DELTA XL AND PICK & GO TECHNOLOGY

- High-acuity bedside monitor that can also go on transport
- Full set of vital parameters, high-acuity parameters via optional Infinity pods
- Integrated Scio gas measurement
- Ergonomic mounting of infusion pumps



ADVANCED WORKSTATION WITH OMEGA MONITORING BRINGING IT-APPLICATIONS TO THE OR

- High-acuity bedside monitor with integrated Pick & Go technology
- 20" touch screen for comprehensive overview of parameters & IT-applications
- Displays patient data from diverse sources: vital signs, ventilation monitoring, lab data, DICOM images & more
- Integrated Scio gas measurement
- Ergonomic mounting of infusion pumps



COMPLETE WORKSTATION WITH DELTA XL MONITORING

- High-acuity bedside monitor that can also go on transport
- Full set of vital parameters, high-acuity parameters via optional Infinity pods
- Integrated Scio gas measurement

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

Locate your Regional Sales
Representative at:
www.draeger.com/contact



**REGION EUROPE CENTRAL
AND EUROPE NORTH**

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

**REGION CENTRAL
AND SOUTH AMERICA**

Dräger Panama Comercial S. de
R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

CANADA

Draeger Medical Canada Inc.
2425 Skymark Avenue, Unit 1
Mississauga, Ontario, L4W 4Y6
Tel +1 905 212 6600
Toll-free +1 866 343 2273
Fax +1 905 212 6601
Canada.support@draeger.com