

LiteStar® Consumables and Accessories

Providing high patient comfort during anesthesia and improving handling for caregivers: The teardrop shape and recessed housing of our anesthesia mask LiteStar support an anatomical fit with effective sealing and provide easier and comfortable grip for your caregiver routines.



Benefits

LiteStar

The Dräger mask LiteStar is intended for single patient use in anesthesia. The new design with the recessed grip supports a perfect control for the physician. In addition, the transparent mask body allows for continuous observation of the patient's condition.

Color-Coding

Every masks size is visualized by a different label color. This labeling supports a quick and easy identification during picking which saves time during the daily routine.

Additional mask accessories

Easy application is supported through additional mask accessories such as a hook ring and headgears for fixation.

Two different headgears in small and large designs offer patients and physicians a soft and elastic support and holding position.

Details



D-37691-2021

Recessed housing side-face




D-37702-2021

Recessed housing front-face



D-44293-2021

Hook ring for LiteStar mask



MT-9654-2006

Headgear, soft rubber



MT-9658-2006

Headgear, soft fabric

Ordering Information

LiteStar anesthesia face mask, disposable, without hook ring	Packaging size	Part number
size 1, small infant	30 pcs.	MP11001
size 2, large infant	30 pcs.	MP11002
size 3, child	30 pcs.	MP11003
size 4, small adult	30 pcs.	MP11004
size 5, medium adult	30 pcs.	MP11005
size 6, large adult	30 pcs.	MP11006
size 7, extra large adult	30 pcs.	MP11007
Hook ring for LiteStar masks		
with 22 mm ID connector (size 3-7)	30 pcs.	MP11008
Headgear, Soft Fabric, disposable		
small	1 pcs.	MP01559
large	1 pcs.	MP01560
Headgear, Soft Rubber, disposable		
small	1 pcs.	MP01561
large	1 pcs.	MP01562

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel: +65 6872 9288
Fax: +65 6259 0398
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06460-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

