



Step up your performance

OXYLOG® 2000 PLUS



Answering the call and stepping up your performance in any situation becomes second nature to you. To those whose lives you touch, you're nothing short of a hero.

From the field to the hospital, you can trust your Oxylog® equipment in any situation

Whether it's responding to an emergency situation or attending to a patient in the Emergency Department, you are asked to step up your performance whilst working under extreme pressure. Skills and experience is what matters. Flexibility is key. Complete confidence in your equipment is essential. The Oxylog 2000 plus can make all the difference.

As the successor of the Oxylog 2000, the Oxylog 2000 plus offers the next step in emergency ventilation. By combining controlled and support modes of ventilation, the Oxylog 2000 plus provides flexibility to adapt the ventilation to your patient's needs. Pressure Support and Non-

Invasive Ventilation are available, enabling support for patients with insufficient breathing and to help prevent intubation at the earliest stage possible. The Oxylog 2000 plus helps you control the details of ventilation and can lead to an improvement in the quality of care.

As a member of the Oxylog family of emergency and transport ventilators, you can rely on the performance of the Oxylog 2000 plus. Proven technology, the result of more than a century of experience in emergency care ventilation, assures consistency in quality and performance. Intuitive operation and an easy-to-read display help you focus on your patient. The Oxylog 2000 plus can become your partner during emergency and transport situations.

OXYLOG® 2000 PLUS

Ventilation modes	VC-CMV, VC-AC, VC-SIMV, SpnCPAP
Options	Pressure Support and Non-Invasive Ventilation
Oxygen Concentration	40%* or 100% (O ₂ Air Mix or 100%)
Tidal Volume	100 mL to 2000 mL
Weight	5.4 kg (including battery) (11.9 lb)
Battery life	4 hours operation time

*patient settings and flows will determine oxygen concentration.

Be more confident...

no matter where the call takes you



Invasive or non-invasive, one solution helps optimize quality of care and patient comfort

A patient's condition can change from one minute to the next. The Oxylog 2000 plus can meet this challenge by putting essential ventilation tools at your fingertips. You can choose the configuration that meets your current treatment needs with the ability to step up that level of performance by adding options designed to improve patient outcomes.

The Oxylog 2000 plus can provide invasive or Non-Invasive Ventilation without having to change your equipment.

Non-Invasive Ventilation as a first line intervention provides a possible alternative in situations where intubation can be avoided; this can lead to a reduced risk of infection and a reduced length of hospital stay.

Control life-threatening situations

The Oxylog 2000 plus is designed to allow direct access to critical settings and patient data via an intuitive user interface. The simple See/Set/Ventilate protocol allows a quick transition from patient assessment to treatment. The large display is divided into three sections; clinical values, airway pressure bar graph and parameter settings, allowing rapid evaluation of the patient's condition during hectic situations.

Proven family platform

Over the past thirty years clinicians have used the Oxylog family of emergency and transport ventilators to provide excellent quality of care to their patients. The Oxylog 2000 plus is the latest addition to the family and is specifically designed for medical professionals who require a device that provides excellent ventilation to patients in emergency situations.



Workload reduction

Ergonomically designed carrying and mounting solutions facilitate quick transitioning from one scene to the next. Robust and reliable power supplies and internal batteries simplify handling and promote flexibility. Standardization of reusable and disposable hose systems facilitate user proficiency and help streamline inventory management.

Maximize device utilization

The Oxylog 2000 plus is backed by comprehensive DrägerService support to enable optimal device performance at all times. Our education and training resources provide support to you and your team to achieve a thorough understanding of the system operation and capabilities.

Step up your performance, step up to the Oxylog 2000 *plus*.

Benefits

- Combining controlled with support ventilation providing the flexibility to adapt to your patient's needs;
- Invasive or Non-Invasive, one solution helps to optimize quality of care and patient comfort;
- The intuitive user interface helps you maintain control of life-threatening situations;
- Proven performance and reliability with the family platform;
- Ergonomically designed transport solutions can reduce your workload;
- Support from DrägerService and also education and training resources maximize device utilization.



MT-4293-2007

Step up your performance

Compact

Specifically designed for use in ambulance services and emergency department.

Access to Information


All parameters available in a high-visibility display enabling confident monitoring, day and night.

Control

Intuitive controls enable operation within seconds of starting the device.

Versatile

Pressure Support and Non-Invasive Ventilation capabilities to address diverse ventilation requirements with one device.



CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

As of August 2015:

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA.