

Dräger SelfTestLung™ Consumables and Accessories

The Dräger SelfTestLung™ combines high-quality materials and an intelligent design to support your clinical workflow. Tested and approved by Dräger for compatibility with a wide range of Dräger ventilators, it can be completely disassembled for easy cleaning. Furthermore, all gas leading parts can be autoclaved at 134°C (273°F).



D-24583-2009

Benefits

Approved for use with a wide range of Dräger ventilators

The Dräger SelfTestLung is designed and tested for use with a wide range of Dräger ventilators for any clinical setting. Its construction includes advanced carbon fiber, silicone, and polysulfone materials which are rugged enough to stand the test of time while delivering the consistent performance you have come to expect from Dräger. This lightweight test lung only weighs in at a mere 190 grams (0.42 pounds).

Complete disassembly for easy cleaning

To make cleaning easier, the Dräger SelfTestLung can be completely disassembled in just one step – a real asset in any busy ICU setting. The carbon-fiber flap can be disinfected to support the hospital's infection control.

Ready to autoclave

With the exception of the carbon fiber flap, the Dräger SelfTestLung can be autoclaved at 134°C (273°F). This ensures the best possible hygiene and minimizes the possibility of crosscontamination. With nosocomial infections on the rise around the world, this feature supports patient safety.

Replacement kit available

Even in the most organized settings, small parts can get lost or go missing. Dräger offers a parts kit that allows you to replace individual missing elements and avoid the cost of replacing the entire lung should any single part become lost or damaged.

Technical Data

Ambient conditions

In storage: Temperature	-10 to +60°C (+14 to +140°F)
In operation: Temperature	+18 to +25°C (+65 to +77°F)
Air pressure	900 to 1,000 hPa (675 to 750 mmHg)
Relative humidity	10 to 65%
Compliance	25 mL/mbar \pm 15% at Vt = 500 mL and PEEP = 0 mbar
Resistance	20 mbar/L/s \pm 15% at 100 L/min
Max. volume	1,000 mL

Ordering Information

Test lung SelfTestLung™, latex-free, reusable, max. 1 L	MP 02 400
Spare parts set SelfTestLung™, with locking brace and attachment bolt	MP 02 401

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Locate your Regional
Sales Representative at:
www.draeger.com/contact

