

Dräger Softbed Mattress Consumables and Accessories

High comfort for our smallest patients: Thanks to the viscoelastic foam core, the biocompatible Dräger Softbed mattress supports patient comfort.



Benefits

Dräger softbed mattress

The skin is the largest human organ. It is of particular importance because it represents a natural barrier to our environment. It is highly sensitive, especially in neonates, and all skin layers may not yet be completely formed in premature infants.

To the point

The design of the Softbed mattress takes into account the special needs of sensitive, immature skin. The soft material minimizes conductive heat losses, which can contribute to the positive development of the infant since cold stress can be avoided.

The inside matters

At the heart of the Softbed is the viscoelastic “memory foam core”. This is a special foam first developed for and used in space travel applications. Its particular characteristic is its ability to be soft and plush but to return to its original shape once the weight has been removed.

High quality standard

The Softbed provides young patients with additional protection: The foam core is manufactured with nontoxic sunflower oil, which poses no health threats. The proven biocompatibility as well as the compliance with the Oeko-Tex standards provide additional, verified safety.

Because time is precious

Handling the Softbed is also a breeze. Dirt is quickly removed from the easy-to-clean cover. If needed, the cover can also be removed entirely and washed. Simply open the contamination protected zipper. Of course, it is possible to also use the Softbed with x-rays.

Ordering Information

Item	Dimensions	Order No.
SoftBed Babyleo ¹	685 x 445 x 30 mm / 27.0 x 17.5 x 1.2 in	2M60390
Softbed Caleo [®] 1	640 x 480 x 40 mm / 25.2 x 18.9 x 1.6 in	MX17012
Softbed Babytherm [®] 1	750 x 470 x 40 mm / 29.5 x 18.5 x 1.6 in	2M21012
SoftBed Isolette [®] C2000/8000/8000+ ¹	740 x 380 x 30 mm / 29.1 x 15.0 x 1.2 in	MP01401
Softbed Resuscitaire w/o scale ^{1, 2}	670 x 540 x 40 mm / 26.4 x 21.3 x 1.6 in	MP01402
Softbed IICS-90 w/o scale ^{1, 2}	770 x 610 x 40 mm / 30.3 x 24.0 x 1.6 in	MP01403
Softbed Resuscitaire with scale ^{1, 2}	650 x 520 x 40 mm / 25.6 x 20.5 x 1.6 in	MP01404
Softbed IICS-90 with scale ^{1, 2}	750 x 590 x 40 mm / 29.5 x 23.2 x 1.6 in	MP01405
SoftBed transport incubator TI500 ¹	635 x 325 x 27 mm / 25.0 x 12.8 x 1.1 in	MP01419

¹Lining: 100% polyurethane, Fabric: 100% polyester

²The foam core is not washable because of the Plexiglas

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION DACH

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
25 International Business Park
#04-20/21 German Centre
Singapore 609916
Tel +65 6308 9400
Fax +65 6308 9401
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377-9100
Fax +507 377-9130
contactcsa@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

